

Comprehensive Mutation Scanning in the Ornithine Transcarbamylase Gene using Hi-Res Melting™ with LCGreen® Plus

J. T. McKinney¹, J. A. Harbour¹, C.R.Ellingson¹, L. Caldovic², M. Lynch², M. Tuchman², S. F. Dobrowolski¹;
¹Idaho Technology, Inc., Salt Lake City, UT, United States, ²Children's National Medical Center, Washington, DC, United States.

Introduction

Ornithine transcarbamylase (OTC, OMIM #311250) deficiency is an X-linked (Xp21.1) inborn error of metabolism involving the urea cycle. OTC deficiency is the most common urea cycle defect (UCD). The algorithm for differential diagnosis of UCD's is complex and labor intensive (Figure 1). Early identification and treatment is critical for a positive prognosis. Pathologic gene mutations in OTC predominantly occur de novo, thus molecular analysis for a "common" mutation is impractical. Herein is presented a rapid method of comprehensively scanning the coding and respective 5' and 3' splice site regions of the OTC gene for mutations.

Methods

Molecular analysis of OTC is useful as an initial reflex to clinical presentation of hyperammonemia. Hi-Res Melting analysis is based on the detection of heteroduplex molecules and combines the use of a new dsDNA binding dye with saturating characteristics (LCGreen Plus) and the LightScanner™ instrument. The dye is included in the PCR reaction and incorporated into newly synthesized DNA during PCR. Following amplification, the samples are transferred directly to the LightScanner for fluorescence monitoring during thermal denaturation (melting). The resulting melting profiles of each sample are then compared to the melting profile of a sequence confirmed wild type sample for the region being queried. Samples with melting profiles that deviate from the profile of the wild type are recovered and sequenced to define the observed variant (Figure 2, panels A-J).

Results

Hi-Res Melting was used to scan the coding and splice regions of the OTC gene in 20 OTC deficient patients. PCR primers were designed to amplify the exons and a minimum of 12 bases of donor/acceptor splice site sequence. A single PCR product was designed to amplify each of the 10 OTC exons in a 96 well microtiter plate. A common amplification protocol was used. Because OTC is X-linked, male samples were co-amplified with a known wild type DNA from a male control, to force heteroduplex formation. Fragment size ranged from 146-266 bp. Nineteen different mutations (14 reported, 5 novel) were identified (Figure 3, novel mutations shown in red).

Conclusion

Hi-Res Melting can be used to rapidly scan for pathologic mutations in the OTC gene. Comprehensive molecular analysis can significantly decrease the time to differential diagnosis of OTC deficiency.

Figure 2. Melting Profiles Observed in 20 OTC Patients




Figure 1. UCD Differential Diagnosis Algorithm





Figure 3. OTC Gene Mutations Identified by Melting



Contact Information

Jason T. McKinney
 jasonm@idahotech.com
 801-736-6354
 www.idahotech.com

Steven F. Dobrowolski
 steven.dobrowolski@idahotech.com
 801-736-6354
 www.idahotech.com